


Edit one analysis through forms, Python or MCP


A Main window: the complete workflow




Detail from A: the two analysis steps

```
▶ *1. Enhanced neuronal signal (root(func=[percentile_normalize]*))
  └─ proc(vars:channel, group:CHANNEL, source:PREVIOUS_STEP) | mat(enabled:True, s
▶ 2. Per-neuron morphology and topology (root(func=neurite_outgrowth_metaxpress))
  └─ proc(vars:channel, group:CHANNEL, source:PREVIOUS_STEP) | nap(enabled:True)
```

B Recorded MCP edits



Connection status in A



C Function controls

Low Percentile: ? 1

High Percentile: ? 99.8

Target Min: ? 0

Target Max: ? 255

D Matching Python code

```
pattern = [
    (percentile_normalize, {
        'low_percentile': 1.0,
        'high_percentile': 99.8,
        'target_min': 0.0,
        'target_max': 255.0,
        'enabled': True
    })
]
```

The same parameters, in the same order, with the same values